



Inhibition of Protein-Protein interaction by a small compound using iMSPR-ProX & COOH-Au chip (2D)

Small compound can inhibit binding of protein 1 to protein 2 by competitively binding to the protein 2 binding site of protein 1. In this application experiment, analysis of inhibition of protein 1 – protein 2 binding was performed at various concentrations of a small compound using iMSPR-ProX.

Materials

- Instrument: iMSPR-ProX
- Sensor chip: COOH-Au chip
- Immobilization Reagent: Amine coupling kit (ACK50)
- Immobilization buffer: Acetate buffer pH5.0 (AB50)
- Regeneration buffer: 10 mM HCl
- Ligand: Protein 2 (6 mg/ml)
- Analyte: Protein 1
- Inhibitor: Small compound

Procedure

Ligand Immobilization

- ① Baseline: 1xHBST is flowed into both channels (ligand channel, reference channel) at a flow rate of 30 μ l/min for

more than 5 minutes to establish a stable baseline.

- ② Inject a 1:1 mixture of NHS/EDC into both channels at a flow rate of 30 μ l/min for 5 minutes followed by 5 minutes of washing in 1xHBST.
- ③ Inject 60 μ g/ml 300 μ l of protein 2 (in Acetate buffer pH5.0) only into the ligand channel at a flow rate of 10 μ l/min for 20 minutes followed by 5 minutes of washing in 1xHBST.
- ④ Inject 200 μ l of Quenching buffer into both channels at a flow rate of 30 μ l/min for 2 minutes followed by 5 minutes of washing in 1xHBST.
- ⑤ Inject 200 μ l of Glycine pH 1.5 into both channels at a flow rate of 30 μ l/min for 2 minutes followed by 30 minutes or more of washing in 1xHBST to stabilize.

Analyte binding

- ① Prepare 200 μ l protein 1 analyte dilutions into the 1xHBST (running buffer) to a concentration of 3.91nM, 7.81nM, 15.6nM, 31.3nM, 62.5nM, 125nM, 250nM, 0 nM.
- ② Introduce 3.91 nM protein 1 at 30 μ l/min for a 3 minutes association time followed by 1xHBST for a 5 minutes dissociation time.
- ③ After binding, the surface was regenerated by injecting 10mM HCl solution at a flow rate of 30 μ l/min for 2 minutes. And then washed with 1xHBST for 6 minutes to stabilize.
- ④ Repeat ②-③ for each concentration

Inhibition




- ① 75 μ l of the small compound was diluted in 1xHBST (150 nM, 1% DMSO) in a series of concentrations (200, 100, 50, 0 μ M) and incubated for more than 30 minutes. At this time, the DMSO concentrations were maintained at 1%.
- ② Introduce mixture of protein 1 (75nM) + small compound 100 μ M at 30 μ l/min for a 3-minute association time followed by 1xHBST for a 5-minute dissociation time.
- ③ After binding, the surface was regenerated by injecting 10mM HCl solution at a flow rate of 10 μ l/min for 3 minutes. And then washed with 1xHBST for 6 minutes to stabilize.
- ④ Repeat ②-③ for each concentration.




Results

R1. Ligand Immobilization




R2. Protein 1-2 binding Curve fitting: one to one binding model



T1. Small compound Kinetic Evaluation

Contents	Value
Immobilization Level	1456.9 RU
B _{max}	14841 RU
K _a (Association rate, 1/M*s)	1.30 x 10 ³
K _d (Dissociation rate, 1/s)	2.49 x 10 ⁻²
K _D (Affinity)	1.91 x10 ⁻⁵ M
Chi ²	29.97

R3. Inhibition



Results summary

The appropriate concentration of protein 1 for inhibition analysis was determined through confirmation of protein 1 and 2 binding and kinetic evaluation.

It was indicated that the higher the concentration of the small compound, the lower the signal of protein 1 binding to Collagen. In particular, it can be observed that most of protein 1 and 2 binding is inhibited at 100 μ M of the small compound.