

A 3D geometric shape composed of several orange rectangular blocks of varying sizes and orientations, creating a complex, angular structure. The text "Surface Plasmon Resonance" is overlaid on the central part of this shape.

**Surface  
Plasmon  
Resonance**

iCLUEB!O

# **Maximultix-ATD on iMSPR-PleX**

Rapid and Precise Discovery  
of Antibody Therapeutics


# Maximultix

**Ready-to-Use,  
Beyond the Dimensions of Analysis  
with Multiplex Pre-array sensor chip**

No more one-by-one verification.

**No complex preparation is required.**

Maximultix is a fully pre-arrayed sensor chip, engineered to be ready the moment you are. Simply load it into our dedicated **iMSPR-Plex** system and transform a single injection into a comprehensive stream of multiplexed data—captured simultaneously and instantly.



iCLUEB!O

# Maximultix-ATD

Antibody Therapeutics Discovery

**The First Innovation.  
Everything for Antibody  
Discovery: Maximultix-ATD.**

Maximultix is purpose-built for therapeutic antibody developers. Now, with just four strategic injections (**1  $\mu$ M, 100 nM / pH7.4, pH6.0**) , you can simultaneously validate structural integrity, efficacy, half-life, and developability. Let the powerful multi-array technology of the Maximultix platform translate your antibody's hidden potential into definitive, high-resolution data.

iCLUEB!O

# 16 Validators

## The 4 Pillars of MMx-ATD Analysis for Antibody Discovery.

### 1. Specificity: The Fundamental Performance.

We precisely capture every moment your antibody identifies and binds to its target, ensuring unwavering specificity.

### 2. Structure & Integrity: A Solid Foundation.

Instantly validate structural integrity and purification status through comprehensive Protein A/G/L and CD64 binding analysis.

### 3. Pharmacokinetics (PK): The Competitive Edge of Longevity.

Predict in-vivo half-life and establish optimal dosing guidelines by analyzing the pH-dependent binding profiles of FcRn.

### 4. Effector Function & Safety: Where Efficacy Meets Security.

Maximize immune activation (ADCC) and minimize side-effect risks through sophisticated CD16 and CD32 binding pattern analysis.



# Analysis Made Easy.

1. Insert the Maximultix-ATD sensor chip into the iMSPR-Plex
2. Open the dedicated ATD App within the Plex OS environment.
3. Inject samples across pH (6.0/7.4) /concentrations (100, 1000 nM)
4. Run the ATD Analysis App for instant, automated data processing and comprehensive reporting.









# Insights Made Clear

## Precision Diagnosis through Automated Multiplexing.

Our platform precisely diagnoses the optimal antibody candidates using automated multiplex-binding curve analysis. From the initial signal to the final report, we go beyond rapid research—integrating multi-dimensional data to deliver the most definitive and comprehensive conclusions.

No.	분류	리간드 (Ligand)	$\Delta$ RU	kon (1/Ms)	koff (1/s)	KD (M)	Status	진단 및 평가 (Diagnosis)
1	Specificity	Target	44.4	6.1E+04	0.00014	2.2E-09	Pass	Criteria: $\Delta$ RU (Target) > $\Delta$ RU (BSA) AND KD < 1E-6 <b>Result: Pass</b> Criteria: Background/NSB signal should be low <b>Result: Background/NSB Check: Clean</b>
2		BSA	7.7	-	-	N.D.	Pass	
3	Structure & Integrity	Protein A	328.9	4.2E+04	0.0004	9.4E-09	Pass	Criteria: (Protein A > 50 or Protein G > 50) AND Protein L > 50 <b>Result: Structural Integrity: Intact (Fc/Fab OK)</b> Criteria: $\Delta$ RU (CD64) > $\Delta$ RU (BSA) <b>Result: Pass</b>
4		Protein G	108.8	4.2E+04	0.00058	1.4E-08	Pass	
5		Protein L	97.2	3.5E+04	0.00071	2E-08	Pass	
6		CD64 (FcγRI)	16.9	4.5E+04	0.0016	3.5E-08	Pass	
7	PK	FcRn (pH 6.0)	5.6	-	-	N.D.	Pass	Criteria: RU(pH6.0) / RU(pH7.4) Ratio > 1 <b>Result: [PK] Excellent (Long Half-life)</b> Criteria: RU < 20 => Clean Profile; RU > 50 => Albumin-Binding Detected <b>Result: Clean Profile</b>
8		FcRn (pH 7.4)	1.4	-	-	N.D.	Pass	
9		HSA (Albumin)	4.8	-	-	N.D.	Pass	
10	Effector	CD16a (V158)	32.1	5.9E+04	0.0054	9.1E-08	Pass	Criteria: $\Delta$ RU (CD16a_V) > $\Delta$ RU (CD16b) AND $\Delta$ RU (CD16a_F) > $\Delta$ RU (CD16b) <b>Result: High ADCC Potential</b> Criteria: CD16b < 20 RU (or < 3×baseline noise) <b>Result: Low Neutrophil Binding (Safe)</b>
11		CD16a (F158)	53.3	4.6E+04	0.02	4.4E-07	Pass	
12		CD16b (NA1)	4.9	-	-	Low Binding	Pass	
13	Safety	CD32a (H)	28.5	-	-	Steady-state	Pass	Criteria: $\Delta$ RU (CD32a H131) > $\Delta$ RU (CD32b) AND $\Delta$ RU (CD32a R131) > $\Delta$ RU (CD32b) <b>Result: Typical IgG1 Binding Profile confirmed.</b>
14		CD32a (R)	26.3	-	-	Steady-state	Pass	
15		CD32b	3.6	-	-	Steady-state	Pass	



Close



# Result guide

## 1. Target Specificity

Criteria	Target > BSA, Target $K_D < 1e^{-6}$ M (Cutoff affinity)
Pass	This indicates that the <b>antibody's variable region (Fv)</b> binds to the <b>target epitope</b> with high affinity. A significantly higher signal compared to the <b>BSA control</b> suggests a low risk of <b>off-target binding</b> , ensuring superior specificity.
Fail	<b>Low binding affinity</b> or excessive <b>non-specific hydrophobic adsorption</b> .

## 2. Structure & Integrity

Criteria	Protein A/G/L > 50RU, CD64 > BSA
Pass	<b>Stable quaternary structure</b> . Strong <b>CD64</b> and <b>Protein A/G/L</b> signals confirm <b>properly folded Fc/Fab</b> and high structural integrity.
Fail	<b>Structural deformation</b> or <b>Fc/Fab degradation</b> . Incomplete folding or aggregation leads to loss of functional binding.

## 3. Pharmacokinetics (PK)

Criteria	FcRn (pH6.0) > FcRn (pH7.4)
Pass	<b>Optimal pH-dependent recycling</b> . High affinity at <b>pH 6.0</b> with rapid dissociation at <b>pH 7.4</b> , confirming superior in-vivo longevity.
Fail	<b>Poor recycling efficiency</b> . Lack of pH-specific binding leads to rapid clearance or lysosomal degradation.

## 4. Effector

Criteria	CD16a(V/F) > CD16b
Pass	<b>Strong CD16a activation</b> . High binding across both <b>V158/F158 polymorphisms</b> ensures potent <b>ADCC activity</b> and <b>broad patient coverage</b> regardless of genetic variation.
Fail	<b>Suboptimal ADCC activation</b> . High <b>fucosylation</b> of the Fc region likely inhibits <b>CD16a binding</b> and immune cell activation.

## 5. Safety

Criteria	CD32a(H/R) > CD32b
Pass	<b>Optimized immune modulation</b> . Effective <b>check-and-balance</b> between activating ( <b>CD32a</b> ) and inhibitory ( <b>CD32b</b> ) receptors ensures controlled immune responses.
Fail	<b>Imbalanced immune response</b> . Excessive <b>CD32b binding</b> may neutralize therapeutic efficacy, while disproportionate <b>activating receptor</b> signals increase the risk of <b>immune toxicity</b> .